



# WAGNER SHAKER

TE-160

Used in the preparation of samples that require a 360° rotational movement.

## Technical Characteristics

### TE-160

- Rotation: 10 to 60 RPM;
- Accuracy:  $\pm 5$  RPM;
- Rotation control: Analog;
- Agitation range: 360°;
- Motor: Reducer 84 Watts;
- Capacity: 08 tests;
- Support: Swivel for 500 ml Stollman bottles, borosilicate;
- Fixing the bottles: Through handles;
- Cabinet: In carbon steel with anti-corrosive treatment and electrostatic painting;
- Dimensions: W = 550 x D = 370 x H = 630 mm;
- Weight: 20 kg;
- Power: 100 Watts;
- Voltage: 220 Volts;
- Accompanies: - 02 extra fuses - Instruction Manual with Warranty Term;
- Note: Does not come with the bottles;

### Benefits and Advantages

- Supports 8 simultaneous samples
- Easy operation, adjusting rotation only
- Fixing the bottles using handles, ensuring process safety
- Carbon steel cabinet with anti-corrosive treatment and electrostatic painting, ensuring longer equipment life and ease of cleaning

## Related Products

